

Playback Switch

SUMMARY

Product Information

Product Name

Playback Switch

Device Category

- | | |
|--|--|
| <input type="checkbox"/> Adapted Toys | <input type="checkbox"/> Kits |
| <input type="checkbox"/> Aids for Daily Living (ADL) | <input type="checkbox"/> LipSyncs |
| <input checked="" type="checkbox"/> Assistive Switches | <input type="checkbox"/> Mounting |
| <input type="checkbox"/> Communication Aids (AAC) | <input type="checkbox"/> Recreation and Leisure |
| <input type="checkbox"/> Computer Access | <input type="checkbox"/> Seating and Positioning |
| <input type="checkbox"/> Environmental Controls | <input type="checkbox"/> Switch Interfaces |
| <input type="checkbox"/> Gaming | <input type="checkbox"/> Writing Aids |
| <input type="checkbox"/> Keyguard | |

User Value Statement

The Playback Switch is an assistive switch that provides audio feedback to the user when activated in addition to functioning as a regular assistive switch.

Designer

Makers Making Change

Device Info

Overview

The Playback Switch is an assistive switch similar to the Interact Switch that also plays a prerecorded message up to 30 seconds long. This device is intended for users of assistive switches who would benefit from or prefer additional audio feedback when using a switch.

Disability Type

Select one or more disability types:

- | | |
|---|-----------------------------------|
| <input checked="" type="checkbox"/> Agility / Dexterity | <input type="checkbox"/> Hearing |
| <input type="checkbox"/> Arthritis | <input type="checkbox"/> Mobility |
| <input type="checkbox"/> Cognitive | <input type="checkbox"/> Mobility |



© 2024 by Makers Making Change

This work is licensed under the CC BY SA 4.0 License: <http://creativecommons.org/licenses/by-sa/4.0>

Playback Switch

SUMMARY

- | | |
|--------------------------------|---------------------------------|
| <input type="checkbox"/> Other | <input type="checkbox"/> SCI |
| <input type="checkbox"/> Pain | <input type="checkbox"/> Vision |

Disability Type Description

The Playback Switch is designed for users of assistive switches who would benefit from or want additional audio feedback on button presses.

How To Use

A secondary user records a message up to 30 seconds long onto the device by holding the black button on the bottom. The Playback Switch can then be connected to a device just like any other assistive switch. The primary user can then press the button, and it will activate the connected device as well as play back the recorded message.

Estimated Cost

The estimated material cost of the device:

- ☒ \$0 - \$10
- ☐ \$11 - \$25
- ☐ \$26 - \$50
- ☐ \$51 - \$100
- ☐ \$101 - \$250
- ☐ \$250+

Attribution

Design by Neil Squire Society/Makers Making Change

Maker Info

Project Skills

- | | |
|--|---|
| <input type="checkbox"/> 3D Printing | <input type="checkbox"/> Other |
| <input type="checkbox"/> Custom PCB | <input type="checkbox"/> Software |
| <input type="checkbox"/> Electronics | <input checked="" type="checkbox"/> Soldering |
| <input type="checkbox"/> Laser Cutting | <input type="checkbox"/> Woodworking |
| <input type="checkbox"/> Mechanics | |



Playback Switch

SUMMARY

Skills Description

This project requires the use of a soldering iron, hot glue gun, and a drill.

Tools Needed

- | | |
|---|--|
| <input type="checkbox"/> 3D Printer | <input type="checkbox"/> Laser Cutter |
| <input checked="" type="checkbox"/> Common Hand Tools | <input checked="" type="checkbox"/> Soldering Iron |
| <input type="checkbox"/> Common Power Tools | <input type="checkbox"/> Specialized Tooling |

Assembly time (hrs)

Estimated Assembly Time: 0.25 hrs

Build Instructions

This project involves opening the commercial button, glueing in an additional button, and soldering in a cable.

Download Link

<https://github.com/makersmakingchange/Playback-Switch/archive/refs/heads/main.zip>

Project Link

<https://github.com/makersmakingchange/Playback-Switch>

License

License

- Everything needed or used to design, make, test, or prepare the Playback Switch is licensed under the [CERN 2.0 Weakly Reciprocal license \(CERN-OHL-W v2\) or later](#).

- Accompanying material such as instruction manuals, videos, and other copyrightable works that are useful but not necessary to design, make, test, or prepare the Playback Switch are published under a [Creative Commons Attribution-ShareAlike 4.0 license \(CC BY-SA 4.0\)](#).

